



# AP Visions


"Lab Connectivity Perfected"

Cost-effective Laboratory Management  
At your fingertips

**Toxicology** This xLab option facilitates complete Toxicology/Pain Management Reporting and Business Rules for both drug screening and drug confirmations, including support for LC/MS instruments and comprehensive reporting.

The Toxicology/Pain Management modules comes complete with:

- ▶ Preconfigured Drug Classes
- ▶ Over 100 Preconfigured common medications
- ▶ Over 100 Preconfigured Metabolites
- ▶ Preconfigured Metabolites to Drug mapping
- ▶ User definable detection window
- ▶ User definable Analyte not detected for reported RX
- ▶ User definable report ordering of Class and Metabolite
- ▶ Comprehensive reporting including
  - Drug Activity Summary Report (Choice of 2)
  - Metabolite Detail report (with and without concentration)
  - Drug Cumulative Report
- ▶ User definable Analyte detected no RX
- ▶ Metabolite active range
- ▶ Specimen Validity

		<b>Laboratory Report</b> Final Report		AP-Visions Laboratory 935 Sundrop Ct Marco Island, FL 34145 Phone: (866)-526-7099 Fax: 706-223-6825 Lab Director: Robert Waechter, Ph.D. CLIA:04C0567999	
<b>PATIENT INFORMATION</b> Patient Name: Louis Pasture MRN: 510488 Patient DOB: 06/17/1938 Sample ID: 100990		<b>REQUESTER INFORMATION</b> Ordering Physician: Donnals Bookman Practice Name: East End Clinic Practice Address: 32145 Seaview Dr Suite 100 Goodland, FL		<b>SAMPLE INFORMATION</b> Specimen type: Urine Collected: 04/25/2019 08:49 AM Received: 04/25/2019 08:49 AM Reported: Ordered: 04/25/2019 08:49 AM	
<b>REPORTED CLINICAL HISTORY - PRESCRIBED MEDICATIONS</b>					
Diazepam [1 Day ago] , Xanax [36 Hours ago]					
<b>CONSISTENT - RESULTS CORRELATE WITH REPORTED CLINICAL HISTORY</b>					
Reported Medication(s)	Anticipated Positive	Amount Detected	Cut-Off na/ml	Outcome	Comments
Diazepam	Diazepam	71	50	POSITIVE	Test result is consistent and expected with reported medication Diazepam.
	Nordiazepam	None Detected	50	NEGATIVE	
	Oxazepam	None Detected	50	NEGATIVE	
	Temazepam	None Detected	50	NEGATIVE	
Xanax	alpha-Hydroxyalprazolam	73	50	POSITIVE	Test result is consistent and expected with reported medication Xanax.
	Alprazolam	91	50	POSITIVE	
<b>INCONSISTENT - ANALYTE DETECTED BUT NO CORRESPONDING MEDICATION REPORTED</b>					
Analyte Detected	Detection Window	Amount Detected ng/mL	Comments		
Paroxetine	Paroxetine is typically detected between 1 Day and 10 Days after last use.	61	Paroxetine was not listed as a prescribed medication.		
<b>REPORTED MEDICATION NOT DETECTED, OR DETECTED OUTSIDE RANGE</b>					
Reported Medication(s)	Metabolite/Parent	Detection Window	Comments		
None	None				
<b>MEDICATIONS REPORTED BUT NOT TESTED</b>					
None					
<b>SPECIMEN VALIDITY</b>					
Test	Outcome	Measured	Reference Range	Comment	
pH	NORMAL	7.10	4.00 to 9.00	None	
Specific Gravity	NORMAL	1.025	1.004 to 1.050	None	
URINE CREATININE	ABNORMAL	19	20 to 1000	None	

This entire laboratory report, including the diagnostic methods used to obtain specific results presented herein, is the property of AP-Visions Laboratories, and is protected by all applicable Federal, State, and Local laws and regulations. Use of this laboratory report and the specific diagnostic results presented herein is limited to the direct licensees of AP-Visions Laboratories, and the improper duplication, publication, and/or dissemination of the format and/or contents of this laboratory report will be considered a violation of Federal, State, and Local laws and regulations. This test was developed and its performance characteristics determined by AP-Visions Laboratories, LLC. It has not been cleared by the U.S. Food and Drug Administration.