

# AP Visions

"Lab Connectivity Perfected"

Cost-effective Laboratory Management  
At your fingertips

**xEMRPLUS** Small laboratories have long needed and wanted an inexpensive way to get their results to their EMR system. For them, purchasing one of the Laboratory Information Systems available today is simply cost prohibitive. **xEmrPlus** is ideal for a variety of small practices (and some large): Pediatric practices, Oncology practices, Urology practices, Waived labs, and any practice that needs an easy to use, cost-effective product.

## What xEmrPlus Is

**xEmrPlus** is a middleware product designed to provide laboratories with the ability to get patient results from their internal lab to their EMR system in a very cost-effective manner.

## What xEmrPlus Is Not AP-Visions never wants to mislead customers.

While **xEmrPlus** has many features common in Laboratory Information Systems, it is not a Laboratory Information System.

## **xEMRPLUS** Software Functionality Includes:

- ▶ A Bidirectional interface for your EMR.  
Send your lab's patient results to your EMR. Elininating manual entry and scanning in your EMR. Receive orders from your EMR (where supported by your EMR).
- ▶ Secure password protected Login  
Protects PHI  
Fully encrypted and secure
- ▶ Age/Gender specific Reference Ranges (including Pediatric)
- ▶ Support for both numeric and alpha results (Like UA's or Flu testing)
- ▶ Manually entered results for your Off-line instruments
- ▶ Support for On-Line (interfaced) instruments
- ▶ Barcode Printing
- ▶ Patient Inquiry to look up historical patient reports and printing patient reports

## OPTIONS

The only option to consider is what online instrument interfaces you need. **xEmrPlus** already comes with everything else!

## WHAT YOU GET

- ▶ **xEmrPlus** software
- ▶ Configuration of your **xEmrPlus** by AP-Visions (Physicians, Locations, Tests, Results, and Ranges)
- ▶ Installation and training
- ▶ One year phone support

**Dashboard - Snapshot of your xEmrPlus**

20 Orders to Match From: eClinical Works Order  
 16 Samples to Review for: CDS Medicor M  
 8 Samples to Review for: Status Plus

**Device List and Status**

Device Name	Device Status
eClinical Works Order	ON LINE
eClinical Works Results	ON LINE
CDS Medicor M	ON LINE
Data Server	ON LINE
Status Plus	ON LINE
Manual Instrument	ON LINE

**Reviewing Numeric Results**

Patient Information: Instrument: CDS Medicor M  
 Select Order to Match: MRN: 50000004 Visit: 30000004 Order: 2015004 Patient: noneyemnonenaxotokebi.radokalewesled DOB: 11/4/2007

Reprint Labels Available for Review 16

Select Run: MRN: 50000004 Lab ID: ??? Specimen: 361 Run D/TL 9/16/2015 8:36:46 AM Instrument: CDS Medicor M

Sample ID: 50000004 Ordering Location: St. Luke's Medical Centre

Patient MRN: 50000004 DOB: 11/04/2007 Sex: M SSN: Mothers Maiden Name: Ordering Physician: Darville, Phyllis

Accepting Tech: User Lab

Analyte	Value	Units	Flag	Reference Range	Comments
<b>Complete Blood Count</b>					
CBC WBC	3.3	x10 <sup>3</sup> /uL	L	9.1 to 34.0	
CBC Lymph	31.5	%		10.0 to 75.0	
CBC Mono	5.2	%		1.0 to 10.0	
CBC Grans	99.4	%	CH	20.0 to 95.0	
CBC Lymph#	3.3	x10 <sup>3</sup> /uL		1.2 to 5.2	
CBC Mono#	0.40	x10 <sup>3</sup> /uL	L	0.40 to 2.00	
CBC Gran#	4.7	x10 <sup>3</sup> /uL		1.5 to 10.0	
CBC RBC	4.18	x10 <sup>6</sup> /uL		4.00 to 6.00	
CBC HGB	17.6	g/dL	H	10.9 to 16.0	
CBC HCT	34.6	%		30.0 to 48.0	
CBC MCV	85	fL		74 to 99	
CBC MCH	28	pg		20 to 40	
CBC MCHC	34	g/dL		30 to 40	
CBC RDW	15.9	%		13.0 to 18.0	
CBC Platelets	273	x10 <sup>3</sup> /uL		150 to 450	
CBC MPV	7.1	fL		6.0 to 9.5	

Clear Screen

Page: 1

**Laboratory Report**  
 Inquiry

Primary Pediatric Care  
 1234 Some Street  
 New York, NY 09912  
 Jonas Salk Dr.  
 Phone: 714-267-8844 Fax: 714-555-1212

TEST	Normal	Abnormal	Flag	Units	Reference Range	Tech
<b>Complete Blood Count</b>						
WBC	31.5	3.3	L	x10 <sup>3</sup> /uL	9.1 to 34.0	
Lymph				%	10.0 to 75.0	
Mono	5.2			%	1.0 to 10.0	
Grans		99.4	CH	%	20.0 to 95.0	
Lymph#	3.3			x10 <sup>3</sup> /uL	1.2 to 5.2	
Mono#		0.40	L	x10 <sup>3</sup> /uL	0.40 to 2.00	
Gran#	4.7			x10 <sup>3</sup> /uL	1.5 to 10.0	
RBC	4.18			x10 <sup>6</sup> /uL	4.00 to 6.00	
HGB		17.6	H	g/dL	10.9 to 16.0	
HCT	34.6			%	30.0 to 48.0	
MCV	85			fL	74 to 99	
MCH	28			pg	20 to 40	
MCHC	34			g/dL	30 to 40	
RDW	15.9			%	13.0 to 18.0	
Platelets	273			x10 <sup>3</sup> /uL	150 to 450	
MPV	7.1			fL	6.0 to 9.5	

**Printable historical Patient Reports**

**Reviewing Alpha Results**

Patient Information: Instrument: Status Plus  
 Select Order to Match: MRN: 50000004 Visit: 30000004 Order: 2015004 Patient: noneyemnonenaxotokebi.radokalewesled DOB: 11/4/2007

Reprint Labels Available for Review 8

Select Run: MRN: 100049 Lab ID: 0006 Specimen: 100049 Run D/TL 1/13/2023 2:47:00 PM Instrument: Status Plus

Sample ID: 100049 Ordering Location: St. Luke's Medical Centre

Patient MRN: 50000004 DOB: 11/04/2007 Sex: M SSN: Mothers Maiden Name: Ordering Physician: Darville, Phyllis

Accepting Tech: User Lab

Analyte	Value	Units	Flag	Reference Range	Comments
<b>UADIP UA DIP</b>					
UADIP Urine Glucose	250	mg/dL		NEGATIVE	
UADIP Bilirubin	Large		A	NEGATIVE	
UADIP Ketone	>=160	mg/dL		NEGATIVE	
UADIP Specific Gravity	1				
UADIP Clarity	Clear			CLEAR	
UADIP pH	>=9.0				
UADIP Protein	30	mg/dL		Negative	
UADIP Urobilinogen	2.0	E.U./dL		0.2 to 1.0	
UADIP Nitrite	Positive		A	Negative	
UADIP Leukocytes	Large		A	Negative	

## CONTACT AP-VISIONS TODAY!

Regardless of the size of your lab or what features you need, AP-Visions has your solution! Contact AP-Visions or your Preferred distributor today for more information and to schedule a demonstration of any of the products AP-Visions has to offer.

**Our sales team stands ready to help!**

### Email Address:

Corporate: Nealf@AP-Visions.Com | Sales Inquires: Sales@AP-Visions.Com  
 Customer Support: CustomerSupport@AP-Visions.Com  
 Reseller Inquiries: Resellers@AP-Visions.Com

### Phone Number: 1-866-526-7099

Corporate: Ext 0 | Sales: Ext 1 | Customer Support: Ext 2 | Resellers: Ext 3